


Copy number variation analysis by next-generation sequencing enhances molecular diagnostic yield of lysosomal storage disorders

Stacey Wong, MS, LCGC
Invitae



Disclosure Information
WORLD*Symposium*TM 2018
Stacey Wong

I am an employee and stockholder of Invitae.

I will not discuss off-label use or investigational use in my presentation.



Introduction

- It is important to confirm the diagnosis of a lysosomal storage disease quickly to initiate management and treatment.
- Enzyme analysis may yield false positive results (i.e., pseudodeficiency alleles) or false negative results, which may not exclude a diagnosis. (PMID: 27293520)

Introduction

- Molecular analysis for diagnosis of lysosomal storage disorders is well established and often used as first-line testing.
 - Read-through sequence analysis alone is insufficient as diagnoses involving CNVs may be missed.
 - Historical methods of copy number variant (CNV) analysis (such as aCGH and MLPA) for LSDs are well established, but use of a second method may delay diagnostic results and treatment. (PMIDs: 22704718, 22848519).
 - Next-generation sequencing (NGS) for CNV detection has not been well characterized in the literature.

Objectives

- Identify molecular causes of disease in individuals with suspected lysosomal storage disorders.
- Assess the contribution of CNVs to positive molecular diagnoses of lysosomal storage disorders.

Methods

- Samples: Genomic DNA samples (n = 444) from blood or saliva
- Next-generation sequencing:
 - Coverage: 350x avg, 50x min
 - + CNV analyses
 - Confirmation of LP/P variants:
 - Sanger or Pacific Biosciences Sequencing
 - Array CGH
- Variant interpretation: ACMG guidelines-based system, SHERLOC (Nykamp *et al.*, *Genet Med.* 2017)


Methods

The Invitae Assay Design

- Coding exons
- Splice sites ± 10 base pairs
- Known pathogenic variants
- Pulldown targets




Copy number analysis (deletion/duplication detection)



Results: Demographics of patients with CNVs identified

- 38 patients identified
- Gender: 20 males, 18 females
- Age range: 0-55 years
- Indications:



Results: Spectrum of CNVs (38/444 patients)

Gene	Classification	Variant details	N	Results outcome
<i>ATP13A2</i>	Pathogenic	Deletion (Exons 17-20)*	1	1 molecular diagnosis
<i>CLN3</i>	Pathogenic	Deletion (Exons 8-9)	9	4 molecular diagnoses, 5 carrier statuses
<i>CLN8</i>	VUS	Duplication (Entire coding sequence)*	1	1 uncertain
<i>CLN8</i>	Pathogenic	Deletion (Entire coding sequence)	1	1 carrier status
<i>CTNS</i>	Pathogenic	Deletion (Exons 1-10)	2	2 molecular diagnoses
<i>DNAJC5</i>	VUS	Duplication (Entire coding sequence)	1	1 uncertain
<i>GAA</i>	Pathogenic	Deletion (Exon 18)	1	carrier status
<i>GALC</i>	Pathogenic	Deletion (Exons 11-17)	8	1 molecular diagnosis, 7 carrier statuses
<i>GLA</i>	Pathogenic	Deletion (Entire coding sequence)	1	1 molecular diagnosis
<i>GLA</i>	VUS	Duplication (Entire coding sequence)*	1	1 uncertain
<i>HEXA</i>	Likely Pathogenic	Deletion (Exons 11-13)*	1	1 molecular diagnosis
<i>IDS</i>	Pathogenic	Deletion (Entire coding sequence)	1	1 molecular diagnosis
<i>KCTD7</i>	VUS	Duplication (Exons 3-4)*	7	7 uncertain
<i>LAMP2</i>	Pathogenic	Deletion (Exon 1)	1	1 molecular diagnosis
<i>LAMP2</i>	VUS	Duplication (Entire coding sequence)*	1	1 uncertain
<i>MFSD8</i>	VUS	Duplication (Exons 2-10)*	1	1 uncertain

* = novel variant

Positive molecular diagnoses: 12 cases with CNVs

Disease	Gene	N	Demographics	CNV	Co-occurring variant
Kufor-Rakeb syndrome	<i>ATP13A2</i>	1	17y male	Deletion (Exons 17-20)	c.1903C>T (p.Q635*)
Neuronal ceroid lipofuscinosis	<i>CLN3</i>	4	5y male, 7y female, 10y female, 13y male	Deletion (Exons 8-9)	c.944dupA (p.H315Qfs*67) c.569dupG (p.A191Sfs*45) 2 homozygotes
Cystinosis	<i>CTNS</i>	2	4y female, 30y female	Deletion (Exons 1-10)	c.206_210delTCCTT (p.I69Rfs*5), 1 homozygote
Krabbe disease	<i>GALC</i>	1	9m male	Deletion (Exons 11-17)	c.1541T>C (p.F514S)
Fabry disease	<i>GLA</i>	1	6y female	Deletion (Entire sequence)	N/A (mosaic Turner syndrome)
Tay-Sachs disease	<i>HEXA</i>	1	1y female	Deletion (Exons 11-13)	homozygote
Hunter syndrome	<i>IDS</i>	1	1y male	Deletion (Entire sequence)	hemizygote
Danon disease	<i>LAMP2</i>	1	16y female	Deletion (Exon 1)	N/A (symptomatic)

Conclusions

- NGS-based CNV analysis is a rapid, efficient, and comprehensive method for diagnosing LSDs.
 - Can shorten turnaround time and reduce cost
- CNV analysis significantly expands the diagnostic yield of LSDs
 - CNVs found in 8.6% (38/444) of patients
 - One-third (12/38) of detected CNVs were clinically actionable
- NGS CNV analysis should be part of first-tier testing for individuals with suspected LSDs.

Acknowledgements



Co-authors: Yuan-Yuan Ho, Laura Murillo, Daniel Beltran, Rachel Harte, Hannah White, Michelle Fox, Tom Winder, Britt Johnson, Invitae Metabolic Team