

INVITAE MICROARRAY PREGNANCY LOSS COLLECTION KIT

IMPORTANT NOTES

Before shipping this kit, the ordering medical professional should complete and include a copy of the pregnancy loss test requisition form (TRF), available for download at www.invitae.com/forms.

1. Label each with the following information:

- Patient's full name, including legal first and last name
- Patient's date of birth (MM/DD/YYYY)
- Patient's sex
- Specimen collection date (today's date)
- Sex of the fetus (if known)

2. Package the specimen(s) appropriately for each type:

For products of conception (POC) fresh tissue

- Place fresh tissue into the empty specimen jar. Pour the contents of the transport media bottle over the tissue in the jar.
- Seal cap tightly to prevent leakage.

For formalin-fixed paraffin-embedded (FFPE) tissue block(s) or unstained slides

- Securely place in the foam insert.

3. Package the kit.

- Place all specimen containers securely in the foam insert.
- Place the foam insert into the biohazard specimen bag with the absorbent pad.
- Place the sealed biohazard bag and TRF in the cardboard box.
- Close the box and place it inside the FedEx Clinical Pak bag.

If you are sending more than one kit in a given day to Invitae, you can include up to three kits in the same FedEx Clinical Pak bag for your convenience.

- Seal the FedEx Clinical Pak bag and apply the included shipping label.
- The kit is now ready for pickup.

4. Schedule a pickup.

- Schedule a pickup with FedEx, either through www.fedex.com or 1-800-GoFedEx (1-800-463-3339).
- If you have questions, please call us at 1-800-436-3037 or email us at clientservices@invitae.com.

STORAGE CONDITIONS

- Store fresh tissue POC sample in a clinical refrigerator (39 °F/4 °C) or at ambient temperature (50–86 °F/10–30 °C).
- Store FFPE blocks and unstained slides at ambient temperature (50–86 °F/10–30 °C).
- Do not freeze. Keep out of direct sunlight.

SHIPPING CONDITIONS

- Ship samples at ambient temperature (50–86 °F/10–30 °C).
- If necessary, ship specimens with a cold pack (i.e., during hot weather conditions).
- Do not freeze.