



To Whom It May Concern:

This letter is to inform those interested that Invitae Corporation's Specimen Collection Kit is composed of a corrugated box containing a specimen collection tube. It is used for transporting a human specimen from a patient back to Invitae Corporation for genetic testing.

The only other components inside are a leak-proof plastic zip lock bag, an absorbent cloth, a piece of paper the patient's physician has filled out, a return shipping label, and a Clinical Lab Pack. The specimen and collection kit are sent via the Clinical Pak envelope back to Invitae Corporation.

Therapak Corporation (address provided below) is the manufacturer, supplier and distributor for the Invitae specimen collection kit. Universal precautions are taken into consideration for human specimen and are included with the specimen collection kit.

Please feel free to contact us with any further questions.

Best regards,

Client Services
Invitae Corporation
1400 16th Street
San Francisco, CA 94103
USA

Therapak Corporation
651 Wharton Drive
Claremont, CA 91711
USA